

兔抗 CNDP2 多克隆抗体

中文名称：兔抗 CNDP2 多克隆抗体

英文名称：Anti-CNDP2 rabbit polyclonal antibody

别名：CNDP dipeptidase 2 (metallopeptidase M20 family) CN2; CPGL; PEPA; HsT2298; HEL-S-13

抗原：CNDP2

储存：冷冻 (-20℃)

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

相关类别：一抗

克隆类型：rabbit polyclonal

技术规格

Background:	CNDP2, also known as tissue carnosinase and peptidase A (EC 3.4.13.18), is a nonspecific dipeptidase rather than a selective carnosinase (Teufel et al., 2003 [PubMed 12473676]).
Applications:	ELISA, IHC
Name of antibody:	CNDP2
Immunogen:	Fusion protein of human CNDP2
Full name:	CNDP dipeptidase 2 (metallopeptidase M20 family)
Synonyms:	CN2; CPGL; PEPA; HsT2298; HEL-S-13
SwissProt:	Q96KP4

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human colorectal cancer
IHC Recommend dilution:	25-100



