

## APOC1 抗原（重组蛋白）

中文名称：APOC1 抗原（重组蛋白）

英文名称：APOC1 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	apolipoprotein C-I
<b>Swissprot:</b>	P02654
<b>Gene Accession:</b>	BC009698
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene is a member of the apolipoprotein C1 family. This gene is expressed primarily in the liver, and it is activated when monocytes differentiate into macrophages. A pseudogene of this gene is located 4 kb downstream in the same orientation, on the same chromosome. This gene is mapped to chromosome 19, where it resides within a apolipoprotein gene cluster. Alternatively spliced transcript variants have been found for this gene, but the biological validity of some variants has not been determined.