

兔抗 ARL5A 多克隆抗体

中文名称：兔抗 ARL5A 多克隆抗体

英文名称：Anti-ARL5A rabbit polyclonal antibody

别名：ARL5; ARFLP5

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ARL5A

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the ARF family of GTP-binding proteins. With its distinctive nuclear/nucleolar localization and interaction with HP1alpha, the protein is developmentally regulated and may play a role (s) in nuclear dynamics and/or signaling cascades during embryonic development. Alternative splicing results in multiple transcript variants encoding different isoforms. This gene has multiple pseudogenes.
Applications:	ELISA, IHC
Name of antibody:	ARL5A
Immunogen:	Synthetic peptide of human ARL5A

Full name:	ADP ribosylation factor like GTPase 5A
Synonyms:	ARL5; ARFLP5
SwissProt:	Q9Y689
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and human tonsil
IHC Recommend dilution:	25-100



