

## 兔抗 ARMC6 多克隆抗体

中文名称：兔抗 ARMC6 多克隆抗体

英文名称： Anti-ARMC6 rabbit polyclonal antibody

别名： R30923\_1

抗原： ARMC6

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	The function of this gene's protein product has not been determined. A related protein in mouse suggests that this protein has a conserved function. Two transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	WB, IF
<b>Name of antibody:</b>	ARMC6
<b>Immunogen:</b>	Synthesized peptide derived from internal of human ARMC6.
<b>Full name:</b>	armadillo repeat containing 6
<b>Synonyms :</b>	R30923_1

<b>SwissProt:</b>	Q6NXE6
<b>WB Predicted band size:</b>	54 kDa
<b>WB Positive control:</b>	Hela cells and COLO cells lysates
<b>WB Recommended dilution:</b>	500-3000
<b>IF positive control:</b>	A549 cells
<b>IF Recommend dilution:</b>	100-500



