

兔抗 CNTN6 多克隆抗体

中文名称：兔抗 CNTN6 多克隆抗体

英文名称： Anti-CNTN6 rabbit polyclonal antibody

别名： NB3

抗原： CNTN6

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The protein encoded by this gene is a member of the immunoglobulin superfamily. It is a glycosylphosphatidylinositol (GPI)-anchored neuronal membrane protein that functions as a cell adhesion molecule. It may play a role in the formation of axon connections in the developing nervous system. Alternative splicing results in multiple transcript variants.
Applications:	WB
Name of antibody:	CNTN6
Immunogen:	Synthesized peptide derived from internal of human CNTN6.

Full name:	contactin 6
Synonyms :	NB3
SwissProt:	Q9UQ52
WB Predicted band size:	114 kDa
WB Positive control:	Jurkat cells and A549 cells lysates
WB Recommended dilution:	500-3000

