

兔抗 ARMC10 多克隆抗体

中文名称：兔抗 ARMC10 多克隆抗体

英文名称： Anti-ARMC10 rabbit polyclonal antibody

别名： SVH; PNAS112; PNAS-112; PSEC0198

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ARMC10

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a protein that contains an armadillo repeat and transmembrane domain. The encoded protein decreases the transcriptional activity of the tumor suppressor protein p53 through direct interaction with the DNA-binding domain of p53, and may play a role in cell growth and survival. Upregulation of this gene may play a role in hepatocellular carcinoma. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 3.
Applications:	ELISA, IHC

Name of antibody:	ARMC10
Immunogen:	Fusion protein of human ARMC10
Full name:	armadillo repeat containing 10
Synonyms :	SVH; PNAS112; PNAS-112; PSEC0198
SwissProt:	Q8N2F6
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	20-100

