

兔抗 ARNTL 多克隆抗体

中文名称：兔抗 ARNTL 多克隆抗体

英文名称： Anti-ARNTL rabbit polyclonal antibody

别名： TIC; JAP3; MOP3; BMAL1; PASD3; BMAL1c; bHLHe5

相关类别： 一抗

抗原： ARNTL

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a basic helix-loop-helix protein that forms a heterodimer with CLOCK. This complex binds an E-box upstream of the PER1 gene, activating this gene and possibly other circadian rhythm-associated genes. Three transcript variants encoding two different isoforms have been found for this gene.
Applications:	ELISA, WB
Name of antibody:	ARNTL
Immunogen:	Fusion protein of human ARNTL
Full name:	aryl hydrocarbon receptor nuclear translocator-like

Synonyms :	TIC; JAP3; MOP3; BMAL1; PASD3; BMAL1c; bHLHe5
SwissProt:	O00327
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	69 kDa
WB Positive control:	Hela cells, mouse heart and stomach tissue
WB Recommended dilution:	200-1000

