

## 兔抗 CNTROB 多克隆抗体

中文名称：兔抗 CNTROB 多克隆抗体

英文名称： Anti-CNTROB rabbit polyclonal antibody

别名： LIP8; PP1221

相关类别： 一抗

抗原： CNTROB

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a centrosomal protein that interacts with BRCA2, and is required for centriole duplication and cytokinesis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	CNTROB
<b>Immunogen:</b>	Fusion protein of human CNTROB
<b>Full name:</b>	centrobin, centrosomal BRCA2 interacting protein
<b>Synonyms :</b>	LIP8; PP1221

<b>SwissProt:</b>	Q8N137
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size:</b>	101 kDa
<b>WB Positive control:</b>	Mouse kidney tissue
<b>WB Recommended dilution:</b>	200-1000

