

兔抗 ARPC1A 多克隆抗体

中文名称: 兔抗 ARPC1A 多克隆抗体

英文名称: Anti-ARPC1A rabbit polyclonal antibody

别 名: ARPC1A; Arc40; SOP2Hs; SOP2L

相关类别: 一抗

储 存: 冷冻(-20℃) 避光

宿 主: Rabbit

抗 原: ARPC1A

反应种属: Human, Mouse

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

This gene encodes one of seven subunits of the huma n Arp2/3 protein complex. This subunit is a member of the SOP2 family of proteins and is most similar to the protein encoded by gene ARPC1B. The similarity bet ween these two proteins suggests that they both may function as p41 subunit of the human Arp2/3 complex that has been implicated in the control of actin polymerization in cells. It is possible that the p41 subunit is involved in assembling and maintaining the structure of the Arp2/3 complex. Multiple versions of the p41 subunit may adapt the functions of the complex to differe



	nt cell types or developmental stages. Alternatively spli ced transcript variants encoding different isoforms have been found for this gene.
Applications:	WB
Name of antibody:	ARPC1A
Immunogen:	Fusion protein of human ARPC1A
Full name:	actin related protein 2/3 complex, subunit 1A
Synonyms:	ARPC1A; Arc40; SOP2Hs; SOP2L
SwissProt:	Q92747
WB Predicted band size:	42 kDa
WB Positive control:	Hela cells
WB Recommended dilution:	500-2000

