

## CBX2 抗原(重组蛋白)

中文名称: CBX2 抗原(重组蛋白)

英文名称: CBX2 Antigen (Recombinant Protein)

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Full length fusion protein

## 技术规格

Full name:	chromobox homolog 2
Synonyms:	M33, CDCA6, SRXY5
Swissprot:	Q14781
Gene Accession:	BC004252
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a component of the polycomb multiprotein complex, which is required to maintain the transcriptionally repressive state of many genes throughout development via chromatin remodeling and modification of histones. Disruption of this gene in mice results in male-to-female gonadal sex reversal. Mutations in this gene are also associated with gonadal dysgenesis in humans. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.