

兔抗 COL5A3 多克隆抗体

中文名称: 兔抗 COL5A3 多克隆抗体

英文名称: Anti-COL5A3 rabbit polyclonal antibody

抗原: COL5A3

储 存: 冷冻(-20℃) 避光

宿 主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:

This gene encodes an alpha chain for one of the low abun dance fibrillar collagens. Fibrillar collagen molecules are tri mers that can be composed of one or more types of alph a chains. Type V collagen is found in tissues containing type I collagen and appears to regulate the assembly of heter otypic fibers composed of both type I and type V collagen. This gene product is closely related to type XI collagen and it is possible that the collagen chains of types V and XI constitute a single collagen type with tissue-specific chain combinations. Mutations in this gene are thought to be responsible for the symptoms of a subset of patients with Ehl ers-Danlos syndrome type III. Messages of several sizes can be detected in northern blots but sequence information ca



	nnot confirm the identity of the shorter messages.
Applications:	IHC, IF
Name of antibody:	COL5A3
Immunogen:	Synthesized peptide derived from internal of human Collag en V $lpha 3$.
Full name:	collagen, type V, alpha 3
Synonyms:	0
SwissProt:	P25940
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100
IF positive control:	Hela cells
IF Recommend dilution:	100-500





