

兔抗 BBS1 多克隆抗体

中文名称: 兔抗 BBS1 多克隆抗体

英文名称: Anti-BBS1 rabbit polyclonal antibody

别名: BBS2L2

抗原: BBS1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Mutations in this gene have been observed in patients with the major form (type 1) of Bardet-Biedl syndrome. The encoded protein may play a role in eye, limb, cardiac and reproductive system development.
Applications:	ELISA, IHC
Name of antibody:	BBS1
Immunogen:	Fusion protein of human BBS1
Full name:	Bardet-Biedl syndrome 1
Synonyms:	BBS2L2
SwissProt:	Q8NFJ9
ELISA Recommended dilution	5000-10000

n:	
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

