

IL20 抗原(重组蛋白)

中文名称: IL20 抗原(重组蛋白)

英文名称: IL20 Antigen (Recombinant Protein)

别 名: IL-20; IL10D; ZCYTO10

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Full length fusion protein

技术规格

Full name:	interleukin 20
Synonyms:	IL-20; IL10D; ZCYTO10
Swissprot:	Q9NYY1
Gene Accession:	BC069311
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a cytokine structurally relate d to interleukin 10 (IL10). This cytokine has been shown to trans duce its signal through signal transducer and activator of transcription 3 (STAT3) in keratinocytes. A specific receptor for this cytok ine is found to be expressed in skin and upregulated dramatically in psoriatic skin, suggesting a role for this protein in epidermal f unction and psoriasis.