

## IL2 抗原(重组蛋白)

中文名称: IL2 抗原(重组蛋白)

英文名称: IL2 Antigen (Recombinant Protein)

别 名: IL-2; TCGF; lymphokine

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Full length fusion protein

## 技术规格

Full name:	interleukin 2
Synonyms:	IL-2; TCGF; lymphokine
Swissprot:	P60568
Gene Accession:	BC066257
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a secreted cytokine that is important for the proliferation of T and B lymphocytes. The receptor of this cytokine is a heterotrimeric protein complex whose gamma chain is also shared by interleukin 4 (IL4) and interleukin 7 (IL7). The expression of this gene in mature thymocytes is monoallelic, which represents an unusual regulatory mode for controlling the precise expression of a single gene. The targeted disruption of a similar gene in



mice leads to ulcerative colitis-like disease, which suggests an essent ial role of this gene in the immune response to antigenic stimuli.