

兔抗 CEP70 多克隆抗体

中文名称：兔抗 CEP70 多克隆抗体

英文名称：Anti-CEP70 rabbit polyclonal antibody

别名：centrosomal protein 70; BITE

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CEP70

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Plays a role in the organization of both preexisting and nascent microtubules in interphase cells. During mitosis, required for the organization and orientation of the mitotic spindle.
Applications:	ELISA, IHC
Name of antibody:	CEP70
Immunogen:	Fusion protein of human CEP70
Full name:	centrosomal protein 70
Synonyms:	BITE
SwissProt:	Q8NHQ1
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human brain
IHC Recommend dilution:	20-100



