

兔抗 CEP72 多克隆抗体

中文名称：兔抗 CEP72 多克隆抗体

英文名称：Anti-CEP72 rabbit polyclonal antibody

别名：centrosomal protein 72

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CEP72

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The product of this gene is a member of the leucine-rich-repeat (LRR) superfamily of proteins. The protein is localized to the centrosome, a non-membraneous organelle that functions as the major microtubule-organizing center in animal cells.
Applications:	ELISA, IHC
Name of antibody:	CEP72
Immunogen:	Fusion protein of human CEP72
Full name:	centrosomal protein 72
SwissProt:	Q9P209
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human thyroid cancer and Human prostate cancer
IHC Recommend dilution:	30-150



