

兔抗 CEP89 多克隆抗体

中文名称: 兔抗 CEP89 多克隆抗体

英文名称: Anti-CEP89 rabbit polyclonal antibody

别名: centrosomal protein 89; CEP123; CCDC123

相关类别: 一抗

抗原: CEP89

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Required for ciliogenesis. Also plays a role in mitochondrial metabolism where it may modulate complex IV activity.
Applications:	ELISA, IHC
Name of antibody:	CEP89
Immunogen:	Synthetic peptide of human CEP89
Full name:	centrosomal protein 89
Synonyms:	CEP123; CCDC123
SwissProt:	Q96ST8
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer

IHC Recommend dilution:

10-50



