

## 兔抗 EIF2S3 多克隆抗体

中文名称：兔抗 EIF2S3 多克隆抗体

英文名称：Anti-EIF2S3 rabbit polyclonal antibody

别名：eukaryotic translation initiation factor 2 subunit gamma; EIF2; EIF2G; MRXSBRK; eIF-2gA; EIF2gamma

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：EIF2S3

反应种属：Human, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is the largest subunit of a heterotrimeric GTP-binding protein involved in the recruitment of methionyl-tRNA(i) to the 40 S ribosomal subunit.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	EIF2S3
<b>Immunogen:</b>	Fusion protein of human EIF2S3
<b>Full name:</b>	eukaryotic translation initiation factor 2 subunit gamma
<b>Synonyms:</b>	EIF2; EIF2G; MRXSBRK; eIF-2gA; EIF2gamma

<b>SwissProt:</b>	P41091
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100



