

兔抗 CEP97 多克隆抗体

中文名称：兔抗 CEP97 多克隆抗体

英文名称： Anti-CEP97 rabbit polyclonal antibody

别名： LRRIQ2; 2810403B08Rik

相关类别： 一抗

抗原： CEP97

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Centrosomal protein of 97 kDa (Cep97), also known as leucine-rich repeat and IQ domain-containing protein 2 (LRRIQ2), is a protein that in humans is encoded by the CEP97 gene. Acts as a key negative regulator of ciliogenesis in collaboration with CCP110 by capping the mother centriole thereby preventing cilia formation. Required for recruitment of CCP110 to the centrosome. Interacts with CALM1, CEP76 and CCP110. Via its interaction with CCP110, may indirectly interact with HERC2 and NEURL4. Interacts with KIF24.
Applications:	ELISA, WB, IHC

Name of antibody:	CEP97
Immunogen:	Fusion protein of human CEP97
Full name:	centrosomal protein 97kDa
Synonyms :	LRRIQ2; 2810403B08Rik
SwissProt:	Q8IW35
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and Human esophagus cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	97 kDa
WB Positive control:	293T, hela, PC3, lovo and A172 cells
WB Recommended dilution:	500-2000



