

## 兔抗 CEP135 多克隆抗体

中文名称：兔抗 CEP135 多克隆抗体

英文名称： Anti-CEP135 rabbit polyclonal antibody

别名： CEP4; MCPH8; KIAA0635

抗原： CEP135

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a centrosomal protein, which acts as a scaffolding protein during early centriole biogenesis, and is also required for centriole-centriole cohesion during interphase. Mutations in this gene are associated with autosomal recessive primary microcephaly-8.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	CEP135
<b>Immunogen:</b>	Synthesized peptide derived from C-terminal of human CEP135.
<b>Full name:</b>	centrosomal protein 135kDa
<b>Synonyms :</b>	CEP4; MCPH8; KIAA0635

<b>SwissProt:</b>	Q66GS9
<b>WB Predicted band size:</b>	133 kDa
<b>WB Positive control:</b>	COLO cells lysate
<b>WB Recommended dilution:</b>	500-3000

