

兔抗 CEP152 多克隆抗体

中文名称：兔抗 CEP152 多克隆抗体

英文名称： Anti-CEP152 rabbit polyclonal antibody

别名： MCPH4; MCPH9; SCKL5

抗原： CEP152

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a protein that is thought to be involved with centrosome function. Mutations in this gene have been associated with primary microcephaly (MCP H4). Alternative splicing results in multiple transcript variants.
Applications:	IHC
Name of antibody:	CEP152
Immunogen:	Synthesized peptide derived from internal of human CEP152.
Full name:	centrosomal protein 152kDa
Synonyms :	MCPH4; MCPH9; SCKL5

SwissProt:	O94986
IHC positive control:	Human lymph noe tissue
IHC Recommend dilution:	50-100

