

兔抗 EIF3J 多克隆抗体

中文名称：兔抗 EIF3J 多克隆抗体

英文名称：Anti-EIF3J rabbit polyclonal antibody

别名：eukaryotic translation initiation factor 3 subunit J; EIF3S1; eIF3-p35; eIF3-alpha

相关类别：一抗

储存：冷冻（-20℃）

抗原：EIF3J

宿主：Rabbit

反应种属：Human, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a core subunit of the eukaryotic initiation factor 3 complex, which participates in the initiation of translation by aiding in the recruitment of protein and mRNA components to the 40S ribosome. There are pseudogenes for this gene on chromosomes 1, 3, and 9. Alternative splicing results in multiple transcript variants encoding different isoforms.
Applications:	ELISA, WB, IHC
Name of antibody:	EIF3J
Immunogen:	Synthetic peptide of human EIF3J
Full name:	eukaryotic translation initiation factor 3 subunit J

Synonyms:	EIF3S1; eIF3-p35; eIF3-alpha
SwissProt:	O75822
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and Human ovarian cancer
IHC Recommend dilution:	30-150
WB Predicted band size:	29 kDa
WB Positive control:	Hela, K562, 293T, Jurkat, 231 and SKOV3 cell lysates
WB Recommended dilution:	200-1000



