

## 兔抗 CERKL 多克隆抗体

中文名称：兔抗 CERKL 多克隆抗体

英文名称：Anti-CERKL rabbit polyclonal antibody

别名：RP26

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CERKL

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein with ceramide kinase-like domains, however, the protein does not phosphorylate ceramide and its target substrate is currently unknown. This protein may be a negative regulator of apoptosis in photoreceptor cells. Mutations in this gene cause a form of retinitis pigmentosa characterized by autosomal recessive cone and rod dystrophy (arCRD).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CERKL
<b>Immunogen:</b>	Synthetic peptide of human CERKL
<b>Full name:</b>	ceramide kinase-like

<b>Synonyms :</b>	RP26
<b>SwissProt:</b>	Q49MI3
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human gastric cancer
<b>IHC Recommend dilution:</b>	50-200

