

兔抗 EIF3K 多克隆抗体

中文名称：兔抗 EIF3K 多克隆抗体

英文名称： Anti-EIF3K rabbit polyclonal antibody

别名： eukaryotic translation initiation factor 3 subunit K; M9; ARG134; PLAC24; PTD001; EIF3S12; HSPC029; MSTP001; PLAC-24; PRO1474; EIF3-

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Rabbit

抗原： EIF3K

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The 700-kD eukaryotic translation initiation factor-3 (eIF3) is the largest eIF and contains at least 12 subunits, including EIF2S12. eIF3 plays an essential role in translation by binding directly to the 40S ribosomal subunit and promoting formation of the 40S preinitiation complex.
Applications:	ELISA, WB, IHC
Name of antibody:	EIF3K
Immunogen:	Fusion protein of human EIF3K
Full name:	eukaryotic translation initiation factor 3 subunit K
Synonyms:	M9; ARG134; PLAC24; PTD001; EIF3S12; HSPC029; MSTP0

	01; PLAC-24; PRO1474; EIF3-p28
SwissProt:	Q9UBQ5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-300
WB Predicted band size:	25 kDa
WB Positive control:	HT29 cell lysate
WB Recommended dilution:	500-2000



