

兔抗 EIF3L 多克隆抗体

中文名称: 兔抗 EIF3L 多克隆抗体

英文名称: Anti-EIF3L rabbit polyclonal antibody

别 名: eukaryotic translation initiation factor 3 subunit L; EIF3EIP; EIF3S11; HSPC021; HSPC025; MSTP005; EIF3S6IP

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Rabbit

抗 原: EIF3L

反应种属: Human, Mouse

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Component of the eukaryotic translation initiation fa ctor 3 (eIF-3) complex, which is required for several steps in the initiation of protein synthesis (PubMed: 17581632, PubMed:25849773, PubMed:27462815). The e eIF-3 complex associates with the 40S ribosome a nd facilitates the recruitment of eIF-1, eIF-1A, eIF-2: GTP:methionyl-tRNAi and eIF-5 to form the 43S pre-initiation complex (43S PIC). The eIF-3 complex sti mulates mRNA recruitment to the 43S PIC and scan ning of the mRNA for AUG recognition. The eIF-3 c



	omplex is also required for disassembly and recyclin g of post-termination ribosomal complexes and sub sequently prevents premature joining of the 40S and 60S ribosomal subunits prior to initiation (PubMed:17581632). The eIF-3 complex specifically targets a nd initiates translation of a subset of mRNAs involved in cell proliferation, including cell cycling, differentiation and apoptosis, and uses different modes of RNA stem-loop binding to exert either translational activation or repression (PubMed:25849773). (Microbial infection) In case of FCV infection, plays a role in the ribosomal termination-reinitiation event leading to the translation of VP2 (PubMed:18056426).
Applications:	ELISA, IHC
Name of antibody:	EIF3L
Immunogen:	Fusion protein of human EIF3L
Full name:	eukaryotic translation initiation factor 3 subunit L
Synonyms:	EIF3EIP; EIF3S11; HSPC021; HSPC025; MSTP005; EIF3 S6IP
SwissProt:	Q9Y262
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	40-200



