

兔抗 CERS1 多克隆抗体

中文名称：兔抗 CERS1 多克隆抗体

英文名称：Anti-CERS1 rabbit polyclonal antibody

别名：LAG1; UOG1; LASS1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CERS1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the bone morphogenetic protein (BMP) family and the TGF-beta superfamily. This group of proteins is characterized by a polybasic proteolytic processing site that is cleaved to produce a mature protein containing seven conserved cysteine residues. Members of this family are regulators of cell growth and differentiation in both embryonic and adult tissues. Studies in yeast suggest that the encoded protein is involved in aging.
Applications:	ELISA, WB, IHC
Name of antibody:	CERS1

Immunogen:	Synthetic peptide of human CERS1
Full name:	ceramide synthase 1
Synonyms :	LAG1; UOG1; LASS1
SwissProt:	P27544
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human ovarian cancer and Human prostate cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	40 kDa
WB Positive control:	Hepg2 and hela cells
WB Recommended dilution:	500-2000



