

## 兔抗 CERT1 多克隆抗体

中文名称：兔抗 CERT1 多克隆抗体

英文名称： Anti-CERT1 rabbit polyclonal antibody

别名： ceramide transporter 1; CERT; GPBP; CERTL; MRD34;

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： CERT1

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a kinase that specifically phosphorylates the N-terminal region of the non-collagenous domain of the alpha 3 chain of type IV collagen, known as the Good pasture antigen. Goodpasture disease is the result of an autoimmune response directed at this antigen. One isoform of this protein is also involved in ceramide intracellular transport. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CERT1

<b>Immunogen:</b>	Fusion protein of human CERT1
<b>Full name:</b>	ceramide transporter 1
<b>Synonyms:</b>	CERT; GPBP; CERTL; MRD34; STARD11; COL4A3BP
<b>SwissProt:</b>	Q9Y5P4
<b>IHC positive control:</b>	Human thyroid cancer and Human gastric cancer
<b>IHC Recommend dilution:</b>	50-200

