

兔抗 EIF4B 多克隆抗体

- 中文名称：兔抗 EIF4B 多克隆抗体
- 英文名称：Anti-EIF4B rabbit polyclonal antibody
- 别名：eukaryotic translation initiation factor 4B; EIF-4B; PRO1843
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：EIF4B
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	Required for the binding of mRNA to ribosomes. Functions in close association with EIF4-F and EIF4-A. Binds near the 5'-terminal cap of mRNA in presence of EIF-4F and ATP. Promotes the ATPase activity and the ATP-dependent RNA unwinding activity of both EIF4-A and EIF4-F.
Applications:	ELISA, IHC
Name of antibody:	EIF4B
Immunogen:	Fusion protein of human EIF4B
Full name:	eukaryotic translation initiation factor 4B
Synonyms:	EIF-4B; PRO1843
SwissProt:	P23588

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human colorectal cancer
IHC Recommend dilution:	100-300



