

兔抗 CETN2 多克隆抗体

中文名称：兔抗 CETN2 多克隆抗体

英文名称：Anti-CETN2 rabbit polyclonal antibody

别名：centrin 2; CALT; CEN2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CETN2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Caltractin belongs to a family of calcium-binding proteins and is a structural component of the centrosome. The high level of conservation from algae to humans and its association with the centrosome suggested that caltractin plays a fundamental role in the structure and function of the microtubule-organizing center, possibly required for the proper duplication and segregation of the centrosome.
Applications:	ELISA, IHC
Name of antibody:	CETN2
Immunogen:	Fusion protein of human CETN2

Full name:	centrin 2
Synonyms:	CALT; CEN2
SwissProt:	P41208
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	200-400

