

兔抗 CFLAR 多克隆抗体

中文名称：兔抗 CFLAR 多克隆抗体

英文名称：Anti-CFLAR rabbit polyclonal antibody

别名：CASH; FLIP; MRIT; CLARP; FLAME; Casper; FLAME1; c-FLIP; FLAME-1; I-FLICE; c-FLIPL; c-FLIPR; c-FLIPS; CASP8AP1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CFLAR

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a regulator of apoptosis and is structurally similar to caspase-8. However, the encoded protein lacks caspase activity and appears to be itself cleaved into two peptides by caspase-8. Several transcript variants encoding different isoforms have been found for this gene, and partial evidence for several more variants exists.
Applications:	ELISA, IHC
Name of antibody:	CFLAR
Immunogen:	Synthetic peptide of human CFLAR

Full name:	CASP8 and FADD-like apoptosis regulator
Synonyms :	CASH; FLIP; MRIT; CLARP; FLAME; Casper; FLAME1; c-FLIP ; FLAME-1; I-FLICE; c-FLIPL; c-FLIPR; c-FLIPS; CASP8AP1
SwissProt:	O15519
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer and Human colon cancer
IHC Recommend dilution:	25-100

