

## 兔抗 CGA 多克隆抗体

中文名称：兔抗 CGA 多克隆抗体

英文名称： Anti-CGA rabbit polyclonal antibody

别名： HCG; LHA; FSHA; GPHa; TSHA; GPHA1; CG-ALPHA

相关类别： 一抗

抗原： CGA

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

The four human glycoprotein hormones chorionic gonadotropin (CG), luteinizing hormone (LH), follicle stimulating hormone (FSH), and thyroid stimulating hormone (TSH) are dimers consisting of alpha and beta subunits that are associated noncovalently. The alpha subunits of these hormones are identical, however, their beta chains are unique and confer biological specificity. The protein encoded by this gene is the alpha subunit and belongs to the glycoprotein hormones alpha chain family. Two transcript variants encoding different isoforms have been found.

	nd for this gene.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	CGA
<b>Immunogen:</b>	Fusion protein of human CGA
<b>Full name:</b>	glycoprotein hormones, alpha polypeptide
<b>Synonyms :</b>	HCG; LHA; FSHA; GPHa; TSHA; GPHA1; CG-ALPHA
<b>SwissProt:</b>	P01215
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human breast cancer and human liver cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	13 kDa
<b>WB Positive control:</b>	HepG2 and 231 cell lysates
<b>WB Recommended dilution:</b>	500-2000



