

兔抗 ELF1 多克隆抗体

中文名称：兔抗 ELF1 多克隆抗体

英文名称：Anti-ELF1 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻 (-20℃)

宿 主：Rabbit

抗 原：ELF1

反应种属：Human, Mouse

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an E26 transformation-specific related transcription factor. The encoded protein is primarily expressed in lymphoid cells and acts as both an enhancer and a repressor to regulate transcription of various genes. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	ELF1
Immunogen:	Fusion protein of human ELF1
Full name:	E74-like factor 1 (ets domain transcription factor)
SwissProt:	P32519
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and human lung cancer

IHC Recommend dilution:

25-100



