

兔抗 ELL 多克隆抗体

中文名称：兔抗 ELL 多克隆抗体

英文名称： Anti-ELL rabbit polyclonal antibody

别名： elongation factor RNA polymerase II

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： ELL

宿主： Rabbit

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Elongation factor component of the super elongation complex (SEC), a complex required to increase the catalytic rate of RNA polymerase II transcription by suppressing transient pausing by the polymerase at multiple sites along the DNA. ELL also plays an early role before its assembly into in the SEC complex by stabilizing RNA polymerase II recruitment/initiation and entry into the pause site. Required to stabilize the pre-initiation complex and early elongation.
Applications:	WB
Name of antibody:	ELL

Immunogen:	Synthesized peptide derived from internal of human ELL.
Full name:	elongation factor RNA polymerase II
Synonyms :	elongation factor RNA polymerase II
SwissProt:	P55199
WB Predicted band size:	68 kDa
WB Positive control:	Jurkat cells and Hela cells lysates
WB Recommended dilution:	500-3000

