

兔抗 ELMO1 多克隆抗体

中文名称：兔抗 ELMO1 多克隆抗体

英文名称：Anti-ELMO1 rabbit polyclonal antibody

别名：CED12; CED-12; ELMO-1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ELMO1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the engulfment and cell motility protein family. These proteins interact with a dedicator of cytokinesis proteins to promote phagocytosis and cell migration. Increased expression of this gene and dedicator of cytokinesis 1 may promote glioma cell invasion, and single nucleotide polymorphisms in this gene may be associated with diabetic nephropathy. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	ELMO1
Immunogen:	Full length fusion protein

Full name:	engulfment and cell motility 1
Synonyms:	CED12; CED-12; ELMO-1
SwissProt:	Q92556
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	30-150

