

兔抗 ELMO2 多克隆抗体

中文名称：兔抗 ELMO2 多克隆抗体

英文名称： Anti-ELMO2 rabbit polyclonal antibody

别名： engulfment and cell motility 2; VMPI; CED12; CED-12; ELMO-2; Ced-12A

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ELMO2

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene interacts with the deicator of cyto-kinesis 1 protein. Similarity to a C. elegans protein suggests that this protein may function in phagocytosis of apoptotic cells and in cell migration. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB, IHC
Name of antibody:	ELMO2
Immunogen:	Synthetic peptide of human ELMO2
Full name:	engulfment and cell motility 2
Synonyms:	VMPI; CED12; CED-12; ELMO-2; Ced-12A
SwissProt:	Q96JJ3

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human thyroid cancer
IHC Recommend dilution:	20-100
WB Predicted band size:	83 kDa
WB Positive control:	Human fetal liver tissue lysate
WB Recommended dilution:	200-1000



