

兔抗 ELOA2 多克隆抗体

中文名称：兔抗 ELOA2 多克隆抗体

英文名称： Anti-ELOA2 rabbit polyclonal antibody

别名： elongin A2; HsT832; TCEB3B; TCEB3L

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ELOA2

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes the transcriptionally active subunit of the SIII (or elongin) transcription elongation factor complex, which also includes two regulatory subunits, elongins B and C. This complex acts to increase the rate of RNA chain elongation by RNA polymerase II by suppressing transient pausing of the polymerase at many sites along the DNA template. Whereas a related protein with similar function, elongin A, is ubiquitously expressed, the encoded protein is specifically expressed in the testis, suggesting it may have a role in spermatogenesis. [provided by RefSeq, Jul 2008]

Applications:	ELISA, IHC
Name of antibody:	ELOA2
Immunogen:	Synthetic peptide of human ELOA2
Full name:	elongin A2
Synonyms:	HsT832; TCEB3B; TCEB3L
SwissProt:	Q8IYF1
IHC positive control:	Human esophagus cancer and Human lung cancer
IHC Recommend dilution:	30-150

