

## 兔抗 ENTPD7 多克隆抗体

中文名称：兔抗 ENTPD7 多克隆抗体

英文名称： Anti-ENTPD7 rabbit polyclonal antibody

别名： LALP1

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： ENTPD7

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

|                          |  |
|--------------------------|--|
| <b>Background:</b>       | LALP1 (lysosomal apyrase-like protein 1), also known as ENTPD7 (ectonucleoside triphosphate diphosphohydrolase 7), is a 604 amino acid multi-pass membrane protein that belongs to the GDA1/CD39 NTPase family. Preferentially hydrolyzing nucleoside 5' -triphosphates, LALP1 and has an order of activity with respect to possible substrates, which is UTP > GTP > CTP. |
| <b>Applications:</b>     | ELISA, IHC   |
| <b>Name of antibody:</b> | ENTPD7   |
| <b>Immunogen:</b>        | Fusion protein of human ENTPD7   |
| <b>Full name:</b>        | ectonucleoside triphosphate diphosphohydrolase 7   |

|                                    |                                       |
|------------------------------------|---------------------------------------|
| <b>Synonyms :</b>                  | LALP1                                 |
| <b>SwissProt:</b>                  | Q9NQZ7                                |
| <b>ELISA Recommended dilution:</b> | 2000-5000                             |
| <b>IHC positive control:</b>       | Human prostate cancer and human brain |
| <b>IHC Recommend dilution:</b>     | 50-200                                |

