

兔抗 EOGT 多克隆抗体

中文名称: 兔抗 EOGT 多克隆抗体

英文名称: Anti-EOGT rabbit polyclonal antibody

别名: EGF domain specific O-linked N-acetylglucosamine transferase; AOS4; AER61; EOGT1; C3orf64

抗原: EOGT

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

This gene encodes an enzyme that acts in the lumen of the endoplasmic reticulum to catalyze the transfer of N-acetylglucosamine to serine or threonine residues of extracellular-targeted proteins. This enzyme modifies proteins containing eukaryotic growth factor (EGF)-like domains, including the Notch receptor, thereby regulating developmental signalling. Mutations in this gene have been observed in individuals with Adams-Oliver syndrome 4. Alternative splicing results in multiple transcript variants.

Applications:	ELISA, IHC
Name of antibody:	EOGT
Immunogen:	Synthetic peptide of human EOGT
Full name:	EGF domain specific O-linked N-acetylglucosamine transferase
Synonyms:	AOS4; AER61; EOGT1; C3orf64
SwissProt:	Q5NDL2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	10-50

