

兔抗 EPB41L2 多克隆抗体

中文名称：兔抗 EPB41L2 多克隆抗体

英文名称：Anti-EPB41L2 rabbit polyclonal antibody

别名：erythrocyte membrane protein band 4.1 like 2; 4.1G; 4.1-G

相关类别：一抗

抗原：EPB41L2

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Required for dynein-dynactin complex and NUMA1 recruitment at the mitotic cell cortex during anaphase (PubMed:23870127).
Applications:	ELISA, IHC
Name of antibody:	EPB41L2
Immunogen:	Synthetic peptide of human EPB41L2
Full name:	erythrocyte membrane protein band 4.1 like 2
Synonyms:	4.1G; 4.1-G
SwissProt:	O43491
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain

IHC Recommend dilution:	50-100
-------------------------	--------

