

兔抗 EPHA10 多克隆抗体

中文名称：兔抗 EPHA10 多克隆抗体

英文名称：Anti-EPHA10 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

抗 原：EPHA10

宿 主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Ephrin receptors, the largest subfamily of receptor tyrosin e kinases (RTKs), and their ephrin ligands are important mediators of cell-cell communication regulating cell attac hment, shape, and mobility in neuronal and epithelial cell s
Applications:	ELISA, IHC
Name of antibody:	EPHA10
Immunogen:	Fusion protein of human EPHA10
Full name:	EPH receptor A10
SwissProt:	Q5JZY3
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human ovarian cancer
IHC Recommend dilution:	25-100

