

兔抗 EPHX1 多克隆抗体

中文名称：兔抗 EPHX1 多克隆抗体

英文名称：Anti-EPHX1 rabbit polyclonal antibody

别名：EPHX

相关类别：一抗

储存：冷冻（-20℃）

抗原：EPHX1

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Epoxide hydrolase is a critical biotransformation enzyme that converts epoxides from the degradation of aromatic compounds to trans-dihydrodiols which can be conjugated and excreted from the body. Epoxide hydrolase functions in both the activation and detoxification of exogenous chemicals such as polycyclic aromatic hydrocarbons. Mutations in this gene cause preeclampsia, epoxide hydrolase deficiency or increased epoxide hydrolase activity.
Applications:	ELISA, IHC

Name of antibody:	EPHX1
Immunogen:	Fusion protein of human EPHX1
Full name:	Epoxide hydrolase 1
Synonyms :	EPHX
SwissProt:	P07099
ELISA Recommended dilution:	500-2000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	10-25

