

兔抗 EPM2A 多克隆抗体

- 中文名称：兔抗 EPM2A 多克隆抗体
- 英文名称：Anti-EPM2A rabbit polyclonal antibody
- 别名：EPM2A, laforin glucan phosphatase; EPM2; MELF
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：EPM2A
- 反应种属：Human, Mouse, Rat
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a dual-specificity phosphatase that associates with polyribosomes. The encoded protein may be involved in the regulation of glycogen metabolism. Mutations in this gene have been associated with myoclonic epilepsy of Lafora. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	EPM2A
Immunogen:	Fusion protein of human EPM2A
Full name:	EPM2A, laforin glucan phosphatase
Synonyms:	EPM2; MELF

SwissProt:	O95278
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human liver cancer
IHC Recommend dilution:	50-300

