

兔抗 EPOR(Ab-426) 多克隆抗体

中文名称：兔抗 EPOR(Ab-426) 多克隆抗体

英文名称： Anti-EPOR(Ab-426) rabbit polyclonal antibody

别名： EPO-R

抗原： EPOR(Ab-426)

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Receptor for erythropoietin. Mediates erythropoietin-induced erythroblast proliferation and differentiation. Upon EPO stimulation, EPOR dimerizes triggering the JAK2/STAT5 signaling cascade. In some cell types, can also activate STAT1 and STAT3. May also activate the LYN tyrosine kinase. Isoform EPOR-T acts as a dominant-negative receptor of EPOR-mediated signaling.
Applications:	WB
Name of antibody:	EPOR(Ab-426)
Immunogen:	Synthesized peptide derived from Internal of human Epo-R.
Full name:	erythropoietin receptor

Synonyms :	EPO-R
SwissProt:	P19235
WB Predicted band size:	55 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000

