

## 兔抗 EPPIN1 多克隆抗体

中文名称：兔抗 EPPIN1 多克隆抗体

英文名称：Anti-EPPIN1 rabbit polyclonal antibody

别名：EPPIN2; EPPIN3; WAP7; WFDC7; dJ461P17.2; CT71; CT72; SPINLW1

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：EPPIN1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

**Background:**

This gene encodes an epididymal protease inhibitor, which contains both kunitz-type and WAP-type four-disulfide core (WFDC) protease inhibitor consensus sequences. Most WFDC genes are localized to chromosome 20q12-q13 in two clusters: centromeric and telomeric. This gene is a member of the WFDC gene family and belongs to the telomeric cluster. The protein can inhibit human sperm motility, and polymorphisms in this gene are associated with male infertility. Read-through transcription also exists between this gene and the downstream WFDC6 (WAP four-disulfide core domain 6) gene.

<b>Applications:</b>	WB
<b>Name of antibody:</b>	EPPIN1
<b>Immunogen:</b>	Fusion protein of human EPPIN1
<b>Full name:</b>	Epididymal protease inhibitor
<b>Synonyms :</b>	EPPIN2; EPPIN3; WAP7; WFDC7; dJ461P17.2; CT71; CT72; SPINLW1
<b>SwissProt:</b>	O95925
<b>WB Predicted band size:</b>	15 kDa
<b>WB Positive control:</b>	Mouse testis tissue
<b>WB Recommended dilution:</b>	500-2000

