

## 兔抗 EPS15 多克隆抗体

中文名称：兔抗 EPS15 多克隆抗体

英文名称：Anti-EPS15 rabbit polyclonal antibody

别名：epidermal growth factor receptor pathway substrate 15; AF1P; AF-1P; MLLT5

相关类别：一抗

储存：冷冻 (-20℃)

抗原：EPS15

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein that is part of the EGFR pathway. The protein is present at clatherin-coated pits and is involved in receptor-mediated endocytosis of EGF. Notably, this gene is rearranged with the HRX/ALL/MLL gene in acute myelogenous leukemias. Alternative splicing results in multiple transcript variants encoding distinct isoforms.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	EPS15
<b>Immunogen:</b>	Synthetic peptide of human EPS15
<b>Full name:</b>	epidermal growth factor receptor pathway substrate 15

<b>Synonyms:</b>	AF1P; AF-1P; MLLT5
<b>SwissProt:</b>	P42566
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer and Human colorectal cancer
<b>IHC Recommend dilution:</b>	50-300



