

## 兔抗 ERAS 多克隆抗体

中文名称： 兔抗 ERAS 多克隆抗体

英文名称： Anti-ERAS rabbit polyclonal antibody

别名： HRAS2; HRASP

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ERAS

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a constitutively active member of the small GTPase Ras protein family. The encoded protein activates the phosphatidylinositol 3-kinase signal transduction pathway in undifferentiated stem cells, but is not expressed in differentiated cells. This gene may be involved in cancer and chemotherapy resistance.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	ERAS
<b>Immunogen:</b>	Synthetic peptide of human ERAS
<b>Full name:</b>	ES cell expressed Ras
<b>Synonyms :</b>	HRAS2; HRASP

<b>SwissProt:</b>	Q7Z444
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human ovarian cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	25 kDa
<b>WB Positive control:</b>	Mouse heart tissue
<b>WB Recommended dilution:</b>	500-2000



