

## 兔抗 ERBB2 (Phospho-Tyr1221/Tyr1222)多克 隆抗体

中文名称: 兔抗 ERBB2 (Phospho-Tyr1221/Tyr1222)多克隆抗体

英文名称: Anti-ERBB2 (Phospho-Tyr1221/Tyr1222) rabbit polyclonal antibody

别 名: NEU; NGL; HER2; TKR1; CD340; HER-2; MLN 19; HER-2/neu

相关类别: 一抗

储 存: 冷冻(-20℃) 避光

宿 主: Rabbit

抗 原: ERBB2 (Phospho-Tyr1221/Tyr1222)

反应种属: Human

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

## 技术规格

cannot bind growth factors. However, it does bind tightly to other ligand-bound EGF receptor family members to form a heterodimer, stabilizing ligand binding and enhancing kinase

-mediated activation of downstream signalling pathways, suc h as those involving mitogen-activated protein kinase and p

This gene encodes a member of the epidermal growth factor (EGF) receptor family of receptor tyrosine kinases. This protein has no ligand binding domain of its own and therefore

hosphatidylinositol-3 kinase. Allelic variations at amino acid



	positions 654 and 655 of isoform a (positions 624 and 625 of isoform b) have been reported, with the most common a llele, Ile654/Ile655, shown here. Amplification and/or overexp ression of this gene has been reported in numerous cancers, including breast and ovarian tumors. Alternative splicing resu lts in several additional transcript variants, some encoding di fferent isoforms and others that have not been fully charact erized.
Applications:	WB, IF
Name of antibody:	ERBB2 (Phospho-Tyr1221/Tyr1222)
Immunogen:	Synthetic peptide of human ERBB2 (Phospho-Tyr1221/Tyr122 2)
Full name:	v-erb-b2 erythroblastic leukemia viral oncogene homolog 2, neuro/glioblastoma derived oncogene homolog (avian) (Phos pho-Tyr1221/Tyr1222)
Synonyms:	NEU; NGL; HER2; TKR1; CD340; HER-2; MLN 19; HER-2/neu
SwissProt:	P04626
WB Predicted band size:	138 kDa
WB Positive control:	MDA cells untreated or treated with EGF
WB Recommended dilution:	500-1000
IF Positive control:	MCF cells
IF Recommended dilution	100-200



