

SOCS6 抗原(重组蛋白)

- 中文名称: SOCS6 抗原(重组蛋白)
- 英文名称: SOCS6 Antigen (Recombinant Protein)

别名: CIS4; SSI4; CIS-4; SOCS4; STAI4; SOCS-4; SOCS-6; STATI4; HSPC060

- 储存: 冷冻(-20℃)
- 相关类别: 抗原

概述

Fusion protein corresponding to a region derived from 384-491 amino acids of human SOCS6

技术规格

Full name:	suppressor of cytokine signaling 6
Synonyms:	CIS4; SSI4; CIS-4; SOCS4; STAI4; SOCS-4; SOCS-6; STATI4; HSPC060
Swissprot:	O14544
Gene Accession:	BC020082
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene contains a SH2 domain and a CIS homolog domain. The protein thus belongs to the cytokine-induced STAT inhibitor (CIS), also known as suppressor of cytokine signaling (SOCS) or STAT-induced STAT inhibitor (SSI), protein family. CIS famil y members are known to be cytokine-inducible negative regulators o f cytokine signaling. The expression of this gene can be induced by GM-CSF and EPO in hematopoietic cells. A high expression level of t



his gene was found in factor-independent chronic myelogenous leuk emia (CML) and erythroleukemia (HEL) cell lines.