

兔抗 ERGIC3 多克隆抗体

中文名称：兔抗 ERGIC3 多克隆抗体

英文名称：Anti-ERGIC3 rabbit polyclonal antibody

别名：ERGIC and golgi 3; Erv46; CGI-54; C2orf47; PRO0989; C20orf47; NY-BR-84; SDBCAG84; dJ477O4.2

相关类别：一抗

储存：冷冻 (-20℃)

宿主：Rabbit

抗原：ERGIC3

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Possible role in transport between endoplasmic reticulum and Golgi.
Applications:	ELISA, WB, IHC
Name of antibody:	ERGIC3
Immunogen:	Fusion protein of human ERGIC3
Full name:	ERGIC and golgi 3
Synonyms:	Erv46; CGI-54; C2orf47; PRO0989; C20orf47; NY-BR-84; SDBCAG84; dJ477O4.2
SwissProt:	Q9Y282
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	50-300
WB Predicted band size:	43 kDa
WB Positive control:	K562 cell lysate
WB Recommended dilution:	500-2000



