

兔抗 ERG 多克隆抗体

- 中文名称: 兔抗 ERG 多克隆抗体
- 英文名称: Anti-ERG rabbit polyclonal antibody
- 别名: p55; erg-3
- 抗原: ERG
- 储存: 冷冻(-20℃)
- 宿 主: Rabbit
- 相关类别: 一抗
- 反应种属: Human, Mouse
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the erythroblast transfor mation-specific (ETS) family of transcriptions factors. All members of this family are key regulators of embryonic development, cell proliferation, differentiation, angiogenes is, inflammation, and apoptosis. The protein encoded by this gene is mainly expressed in the nucleus. It contains an ETS DNA-binding domain and a PNT (pointed) domai n which is implicated in the self-association of chimeric o ncoproteins. This protein is required for platelet adhesion to the subendothelium, inducing vascular cell remodeling. It also regulates hematopoesis, and the differentiation an
-------------	--



	d maturation of megakaryocytic cells. This gene is involve d in chromosomal translocations, resulting in different fus ion gene products, such as TMPSSR2-ERG and NDRG1-ER G in prostate cancer, EWS-ERG in Ewing's sarcoma and F US-ERG in acute myeloid leukemia. Multiple alternatively spliced transcript variants encoding different isoforms hav e been identified.
Applications:	ELISA, WB, IHC
Name of antibody:	ERG
Immunogen:	Fusion protein of human ERG
Full name:	V-ets erythroblastosis virus E26 oncogene homolog (avian)
Synonyms :	p55; erg-3
SwissProt:	P11308
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colon cancer and Human lung cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	54 kDa
WB Positive control:	293T cells
WB Recommended dilution:	500-2000







